

# Supernal Sleep System

All models available in Twin 80, Full 80, Queen & Dual King.



Features	Supernal 1 Matt Mover	Supernal 2 Basic	Supernal 2 Recliner Plus	Supernal 3 Basic	Supernal 3 Hi-Low	Supernal 5
Head & Foot	X	✓	✓	✓	✓	✓
Hi-Low	✓	X	X	✓	✓	✓
Trendelenburg	X	X	X	X	X	✓
Reverse Trendelenburg	X	X	X	X	X	✓
Cardiac Chair	X	X	X	X	X	✓
Wallhugging	X	X	✓	X	✓	X
European Style Head Tilt	X	X	✓	X	✓	X
Twin Capacity	400 lbs	400 lbs	400 lbs	400 lbs	400 lbs	300 lbs
Full/Queen Capacity	500 lbs	500 lbs	500 lbs	500 lbs	500 lbs	300 lbs
Standard Hand Control	Wireless	Wireless	Wireless	Wireless	Wireless	Wireless
Massage	X	X	✓	X	✓	✓
Hoyer Compatible*	w/ Casters**	Low profile	Low profile	w/ Casters**	w/ Casters**	Low profile
Environmental Control Unit (ECU) Support	X	X	X	X	X	X

\*\* 5" Locking Casters enables use of low profile Hoyer

Technical Specs	Supernal 1 Matt Mover	Supernal 2 Basic	Supernal 2 Recliner Plus	Supernal 3 Basic	Supernal 3 Hi-Low	Supernal 5
Deck Height	7"	14"	11.5"	10"	10.5"	12"
Range	10"	N/A	N/A	10"	10"	12"
Head Angle	N/A	60°	65°	65°	65°	60°
Foot Angle	N/A	45°	35°	35°	35°	35°
Trendelenburg Angle	N/A	N/A	N/A	N/A	N/A	10°
Reverse Trendelenburg Angle	N/A	N/A	N/A	N/A	N/A	10°
Cardiac Chair Angle	N/A	N/A	N/A	N/A	N/A	70°
Twin 80 Frame Weight*	168 lbs	95 lbs	151 lbs	215 lbs	263 lbs	186 lbs
Full 80 Frame Weight*	235 lbs	120 lbs	199 lbs	256 lbs	335 lbs	232 lbs
Queen Frame Weight*	247 lbs	130 lbs	204 lbs	266 lbs	340 lbs	235 lbs
Current	DC	DC	DC	DC	DC	DC
Duty Cycle	10%	10%	10%	10%	10%	10%

\*Uncrated, without mattress

v2 - 11.22.21