

Rehab Line Standard Beds



| Bed Features | Night Rider | New Valiant | Companion |
|--|-------------|-------------|-------------|
| Deck Height | 12" | 12" | 12" |
| Range | 12" | 12" | N/A |
| Head & Foot | ✓ | ✓ | ✓ |
| Hi-Low | ✓ | ✓ | ✗ |
| Trendelenburg | ✗ | ✓ | ✗ |
| Reverse Trendelenburg | ✗ | ✓ | ✗ |
| Cardiac Chair | ✗ | ✓ | ✗ |
| Twin 80 & 84 | ✓ | ✓ | ✓ |
| Full 80 & 84 | ✓ | ✓ | ✓ |
| Dual King | ✓ | ✓ | ✓ |
| Special Size 48x80 | ✓ | ✓ | ✓ |
| Capacity | 300 lbs | 300 lbs | 300 lbs |
| Standard Hand Control | Wired | Wired | Wired |
| Deluxe Upgrade Available | ✓ | ✓ | ✓ |
| Hoyer Compatible* | Low profile | Low profile | Low profile |
| Environmental Control Unit (ECU) Support | ✓ | ✓ | ✓ |

| Technical Specs | Night Rider | New Valiant | Companion |
|-----------------------------|-------------|-------------|-----------|
| Head Angle | 60° | 60° | 60° |
| Foot Angle | 35° | 35° | 35° |
| Trendelenburg Angle | N/A | 10° | N/A |
| Reverse Trendelenburg Angle | N/A | 10° | N/A |
| Cardiac Chair Angle | N/A | 70° | N/A |
| Twin 80 & 84 Frame Weight** | 216 lbs | 216 lbs | 168 lbs |
| Full 80 & 84 Frame Weight** | 270 lbs | 270 lbs | 213 lbs |
| 48x80 Frame Weight** | 270 lbs | 270 lbs | 213 lbs |
| Queen Frame Weight** | 279 lbs | 279 lbs | 222 lbs |
| Current | AC | AC | AC |
| Duty Cycle | 20% | 20% | 20% |

** Crated, without mattress

v2 - 11.22.21