

# Rehab Line Standard Beds



Bed Features	Night Rider	New Valiant	Companion
Deck Height	12"	12"	12"
Range	12"	12"	N/A
Head & Foot	✓	✓	✓
Hi-Low	✓	✓	✗
Trendelenburg	✗	✓	✗
Reverse Trendelenburg	✗	✓	✗
Cardiac Chair	✗	✓	✗
Twin 80 & 84	✓	✓	✓
Full 80 & 84	✓	✓	✓
Dual King	✓	✓	✓
Special Size 48x80	✓	✓	✓
Capacity	300 lbs	300 lbs	300 lbs
Standard Hand Control	Wired	Wired	Wired
Deluxe Upgrade Available	✓	✓	✓
Hoyer Compatible*	Low profile	Low profile	Low profile
Environmental Control Unit (ECU) Support	✓	✓	✓

Technical Specs	Night Rider	New Valiant	Companion
Head Angle	60°	60°	60°
Foot Angle	35°	35°	35°
Trendelenburg Angle	N/A	10°	N/A
Reverse Trendelenburg Angle	N/A	10°	N/A
Cardiac Chair Angle	N/A	70°	N/A
Twin 80 & 84 Frame Weight**	216 lbs	216 lbs	168 lbs
Full 80 & 84 Frame Weight**	270 lbs	270 lbs	213 lbs
48x80 Frame Weight**	270 lbs	270 lbs	213 lbs
Queen Frame Weight**	279 lbs	279 lbs	222 lbs
Current	AC	AC	AC
Duty Cycle	20%	20%	20%

\*\* Crated, without mattress

v2 - 11.22.21